



# Chemical Transfer Unit

## Single product




Transfer product up to 5 gpm (18.9 lpm)

Available with Santoprene, Viton, or Kalrez pump seals

Operate unit by using manual air button, or timer and solenoid valve

Includes suction and discharge hoses



**Options:** Add option codes to item number as shown below

	Pump Seal Material	Operating Type	
CT	Santoprene (1)	Manual Air Button (A)	5
	Viton (2)	Timer and Solenoid Valve (D)	
	Kalrez (5)		

Examples:

- CT1A5 (Single product transfer unit with Santoprene pump seals and manual air button)
- CT2D5 (Single product transfer unit with Viton pump seals and timer and solenoid valve)
- CT5A5 (Single product transfer unit with Kalrez pump seals and manual air button)