



# Chemical Transfer Unit

## Single product




Transfer product up to 3.5 gpm (13.2 lpm)

Available with Santoprene, Viton, or Kalrez pump seals

Operate unit by using manual air button, or timer and solenoid valve

Includes suction and discharge hoses



**Options:** Add option codes to item number as shown below

	Pump Seal Material	Operating Type	
CT	Santoprene (1)	Manual Air Button (A)	5
	Viton (2)	Timer and Solenoid Valve (D)	
	Kalrez (5)		

Examples:

- CT1A5 (Single product transfer unit with Santoprene pump seals and manual air button)
- CT2D5 (Single product transfer unit with Viton pump seals and timer and solenoid valve)
- CT5A5 (Single product transfer unit with Kalrez pump seals and manual air button)